


## TREATMENT FOR BRAIN FATIGUE

Effective points for applying moxa:

- **MS1 Ezhong**
- **Ex-HN3 Yintang**
- **GB14 Yangbai**
- **Ex-HN5 Taiyang**

### ATTENTION

- Turn down/off air conditioner to stabilize air flow.
- Take care not to startle patient with the sound of lighter.
- Make sure moxa is adhered securely when you put it laterally on the face.
- Slowly twist moxa when you detach from skin surface.




## TREATMENT FOR EYESTRAIN

Effective points for applying moxa:

- **Ex-HN4 Yuyao**
- **GB14 Yangbai**
- **ST2 Sibai**
- **Ex-HN5 Taiyang**

### ATTENTION

- Turn down/off air conditioner to stabilize air flow.
- Take care not to startle patient with the sound of lighter.
- Use extra caution when adhering to the eyebrow.
- Slowly twist moxa when you detach from skin surface.



## TREATMENT FOR TOOTHACHE

Effective points for applying moxa:

- **GB3 Shangguan**
- **ST6 Jiache**
- **ST2 Sibai**
- **ST5 Daying**

### ATTENTION

- Turn down/off air conditioner to stabilize air flow.
- Take care not to startle patient with the sound of lighter.
- Make sure moxa is adhered securely when you put it laterally on the face.
- Slowly twist moxa when you detach from skin surface.
- Apply moxa below nose after smoke reduces to avoid breathing in.

